

BraceID[®]

Comfort

Let's have a look!





Torso



■ Universal Rigid Neck Orthosis	4
■ Cervical Collars	5
■ Lumbo Sacral Orthosis	6
■ Contoured Back Support	7
■ Double Cross Back Support	7
■ SpineBoard: Dorsal Splint with Back Support	8
■ Abdominal Bandage	10

Upper Extremity



■ Universal Comfort Arm Sling	12
■ OrthoSling	12
■ Comfort Shoulder Immobiliser	13
■ Shoulder ABDuction Pillow 15°	14
■ ROM Elbow Orthosis	15
■ Elbow Strap	16
■ D-Ring Wrist Brace & LONG D-Ring Wrist Brace	17
■ D-Ring Wrist-Thumb Brace	18
■ Thumb Spica Splint	19
■ Stack Finger Splints	19

Lower Extremity



■ Functional 4-Point Knee Orthosis	21
■ Universal Post-Op ROM Knee Orthosis	22
■ COOL ROM Knee Orthosis	23
■ Tri-Panel Knee Immobilizer	24
■ Knee Support with Hinge	25
■ Knee Support with Stays	26
■ Post Trauma Knee Dressing (PTKD)	27
■ Tri-Step Ankle Lace-Up Brace	28
■ Universal Ankle Stirrup Brace	29
■ Ankle Wrap	30
■ Shield Walkers Hi Top	31
■ LOW Shield Walker	32
■ Walker Hi Regular	33
■ Short Regular Walker	34
■ ROM Walkers (Standard liner/AIR liner)	35
■ KiDS Walker	36

TORSO



Comfort fit

Universal Rigid Neck Orthosis



Features and Benefits

- Universal Collar, easy to adjust in height and circumference.
 - Height adjustment with (green) button system at lateral, front side.
 - Circumference adjustment with hook and loop closure.
- Pivoting chin support.
- High quality plastic with flexible edges.
- Terry towel foam liner provides extra cushioning comfort.
- Universal size that covers a large part of the patient population.

Indications

- Post-operative immobilisation
- Traumatic neck injuries with significant instabilities
- Severe cervical syndrome
- Cervical herniated disc
- Stable fractures of the cervical vertebrae



Neck circumference

25 - 48 cm

Description

Rigid Neck Orthosis

Code No.

12-0082-1

Product Composition

Polyester, Nylon, PU foam, Plastic



Features and Benefits

- Pre-shaped, contoured foam for a comfortable, anatomical fit.
- High back side provides superior support in extension.
- Rounded foam edge, facing towards the skin, offers extra wearing comfort.
- Additional washable, cotton stockinet cover provides improved patient hygiene.
- Variable hook and loop closures suitable for a wide range of neck circumferences.
- Available in a selection of heights to suit a wide range of patients.

Indications

- Acute or chronic painful cervical syndrome
- Whiplash accidents
- Rheumatoid Arthritis



Neck circumference	Neck height	Code No.
22 - 48 cm	6 cm	11-0086-6
28 - 49 cm	8 cm	11-0086-8
28 - 52 cm	10 cm	11-0086-10
28 - 52 cm	12 cm	11-0086-12

Product Composition

PU foam, Cotton, Nylon, Lycra



Features and Benefits

- Modular system with removable posterior panel and removable side pads.
 - Allows the practitioner complete customisation for a treatment according severity of the patient's problem or during progress in rehab.
 - Provides comfort and patient compliance.
- The BraceID® Lumbo Sacral Orthosis, with rigid anterior and posterior moulded frames, extends from the sacrococcygeal junction to the T9-Vertebrae.
 - Ergonomic posterior panel for more targeted motion control in the area S1-T9.
 - Lateral strength is provided by rigid lateral pads. These pads produce intra-cavitary pressure to reduce load on intervertebral discs. The compression and lateral strength is created by a mechanical pulley system that needs to be adjusted by the patient (pulling or releasing the cords).
- The concept allows patients to actively participate in the treatment protocol.
- By increasing the intra-cavitary pressure, the patient can reduce the load (pressure) on the spine.
- Size: cut and adjust - S to 5XL.

Indications

- Lumbar mechanical back pain
- Spinal stenosis
- Back/spine trauma
- Spondylolisthesis
- Facet syndrome
- Post-Op stabilisation protocol, following Laminotomy, Foraminotomy, Anterior Laproscopic fusion and Multilevel Decompression

Waist circumference	Size	Code No.
70 - 170 cm	Universal	12-0081-3

Product Composition

Polyester, Polyamide, Spandex, PU Foam, Neoprene and Plastic



Easy pull system



Stabilising back support



Contoured Back Support

Features and Benefits

- Contoured elastic support with anatomical fit.
- Extra side pulls at the back emphasise the lordosis and provide optimal support for the back.
- Four semi-rigid aluminium stays in the back part with lordosis contour.
- Two spring stays at the front for extra support.
- Hook and loop closure at the front with an additional finger loop for easy application.
- Front height 18 cm / 26 cm at the back.



Waist circumference	Size	Code No.
67 - 80 cm	Small	11-0081-S
78 - 92 cm	Medium	11-0081-M
89 - 105 cm	Large	11-0081-L
103 - 130 cm	X-Large	11-0081-XL

Indications

- Low back pain
- Sprains and strains
- Overuse injuries

Product Composition

Polyester, Nylon, Spandex, PU Foam, Aluminium and Plastic



Double Cross Back Support

Features and Benefits

- High quality rib elastic for an improved skin aeration.
- Dorsal double cross strap for an optimal back support.
- Four semi-rigid aluminium stays in the back part with lordosis contour.
- Two plastic stays at the front for extra support.
- Hook and loop closure at the front with an additional finger loop for easy application.
- Front height 19 cm / 26 cm at the back.

Indications

- Low back pain
- Sprains and strains
- Overuse injuries



Waist circumference	Size	Code No.
69 - 80 cm	Small	11-0082-S
80 - 96 cm	Medium	11-0082-M
96 - 112 cm	Large	11-0082-L
112 - 130 cm	X-Large	11-0082-XL
130 - 148 cm	XX-Large	11-0082-XXL

Product Composition

Polyester, Nylon, Spandex, PU foam, Aluminium, Plastic

SpineBoard: Dorsal Splint with Back Support



The SpineBoard orthosis consists of a padded plastic splint reinforced with two aluminum stays and inserted into a textile lumbar support bandage. Additional shoulder straps pull the patient's body backwards, against the splint.

Wearing the SpineBoard activates the back muscles. The upper part of the body straightens itself, counteracting the hunched back that is typical to osteoporosis. This splint should be individually adjusted by an orthotist when worn the first time. The BraceID® SpineBoard has been designed to enable the patient to breathe normally and to move arms and shoulders without restriction. The device has been adapted specifically for senior citizens, with limited mobility of the limbs, as a target group.



Features and Benefits

- The plastic splint's contoured S-shape corresponds to the normal lumbar lordosis and thoracic kyphosis.
 - A vertical groove in the centre alleviates pressure on the bony protrusions of the processus spinosus.
 - Two thin, shaped aluminium bars at each side of this groove provide support to the plastic splint. The splint has been fitted with extra-soft padding strips inside for additional comfort.
 - The splint is covered in fabric so that it can be worn discretely.

- Each shoulder strap has a metal buckle with a roller system to make application easier for patients in all categories. Moreover, the ends of each strap with the removable Y-tab closure system taper outwards to prevent them from accidentally sliding through the rings.⁽¹⁾
- The back support has also extra features to help elderly persons with limited dexterity.
 - The back bandage has wide finger pockets to facilitate application.
 - Extra finger loops have been attached at each end of the side pull straps to enable patients to grasp them more easily at the side of the body.
- These elastic side pull straps are applied in a V-shape at the back of the support bandage to keep the spinal splint as close to the body as possible.
- Less SKU's results in a reduced stock level. Combine splint and waist bandage according the patient's body size and shape:
 - 4 lengths for Spinal Splints (height: 40-46-51-56 cm)
 - 5 Back supports (waist circumference: from 50 to 140 cm)

Indications



- Osteoporotic vertebral body fractures of the thoracic and/or lumbar spine
- Juvenile osteochondrosis
- Kyphosis with chronic backache
- Osteoporosis



(1)



Add A to B

	Back Support		1x A + 1x B	Dorsal Splint		
	A	Size		Size	B	
50 - 64 cm	12-0084-XS	X-Small		X-Small	12-0083-40	40 cm
63 - 75 cm	12-0084-S	Small		Small	12-0083-46	46 cm
73 - 90 cm	12-0084-M	Medium		Medium	12-0083-51	51 cm
88 - 105 cm	12-0084-L	Large		Large	12-0083-56	56 cm
103 - 125 cm	12-0084-XL	X-Large				
120 - 140 cm	12-0084-XXL	XX-Large				

Product Composition : SpineBoard
Polyester, Nylon, Spandex and PU foam

Abdominal Bandage



Features and Benefits

- A lightweight and breathable rib elastic increases wearing comfort.
- Hook and loop closing system facilitates application and removal.
- Elastic fabric is also hook receptive, which makes size choice easier.
- Washable for good patient comfort and hygiene.
- Latex free.

Indications


- Post-surgical interventions
- Support after cosmetic abdominal surgery or liposuction

Waist circumference	Size	Code No.
65 - 75 cm	X-Small	11-0085-1
70 - 86 cm	Small	11-0085-2
80 - 95 cm	Medium	11-0085-3
90 - 105 cm	Large	11-0085-4
100 - 125 cm	X-Large	11-0085-5
120 - 140 cm	XX-Large	11-0085-6

Product Composition
Polyester, Nylon, PU foam, Spandex



UPPER
EXTREMITY



Hmmm...
Enjoy your rehab!

Universal Comfort Arm Sling



Features and Benefits

- BraceID® Comfort Arm Sling is a lightweight, extremely soft and comfortable arm support.
- The stretch fabric of the envelope adjusts itself easily to the contour and size of the arm.
- Pockets inside the envelope allow to store personal items and give easy access to it for the one handed person.
- A wide and soft neck strap provides extra comfort for neck and shoulder.
- Closure system with Y-tab hook (strap can be easily trimmed to correct length).
- Universal size, fits left or right.

Indications

- Use following trauma injuries, and chronic or neurological conditions.
- Supports the arm with elbow in bended position

Forearm length	Size	Code No.
35 - 45 cm	Universal	11-0064-U

Product Composition

Polyester, Nylon, PU Foam, Metal

OrthoSling



The OrthoSling is an easy and quickly to apply immobilisation system for shoulder and elbow.

Features and Benefits

- Soft, non-stretch and comfortable PU foam strap.
- Easy to handle Hook Closure system, fixed at one side, and removable (y-tab) at the other side. This allows a quick length adjustment of the strap.
- 190 cm length
- Hospital pack of 5 units per box.
- Versatile in application

Indications

- Post-Op or Post-Trauma immobilisation of the upper extremity.
- Support of the arm after a cast application.



Type	Strap length	Code No.
Five Unit Box	190 cm each	11-0061-5

Product Composition

Polyester, Nylon, PU Foam



Comfort Shoulder Immobiliser

Features and Benefits

- Extra soft neck pad provides additional patient comfort.
- In length adjustable neck and waist straps with a Y-tab hook system facilitates secure closure and correct clinical positioning.
- A press button closure system to open the envelope allows early elbow mobilisation and personal hygiene without complete sling removal.
- The particular design of waist strap attachment on the sling enables wrist movement at all times.
- Soft thumb loop effectively maintains sling placement.
- **COMFORT PLUS:** extra convenient pockets on the sling provide benefits to the one-handed patient.



Indications

- Post-trauma immobilisation
- Post-Surgical immobilisation
- Can be used following shoulder surgery, fracture, and or surrounding soft tissue injury

Forearm length	Size	Comfort	Comfort Plus
14 - 19 cm	X-Small	11-0065-XS	
19 - 25 cm	Small	11-0065-S	11-0068-S
25 - 30 cm	Medium	11-0065-M	11-0068-M
30 - 35 cm	Large	11-0065-L	11-0068-L
> 35 cm	X-Large	11-0065-XL	11-0068-XL

Product Composition

Nylon, Polyester, PU foam, Plastic, Metal



Shoulder ABDuction Pillow 15°



Features and Benefits

- The three-point closure system avoids pressure in the neck and is easy to apply. In addition, extra soft neck and arm pit padding provide extra comfort.
- Click buckles ensure safe fixation and correct clinical positioning of the arm.
- Strap length can be easily adjusted by trimming them and repositioning Y-tab hook straps.
- The soft envelope can be opened quickly for arm mobilization or personal hygiene without having to completely remove sling or pillow.
- Instructions for left or right side application are simple due to the colour code of neck and arm pit straps. (green/black)
- One extra hygiene (red) replacement strap in each kit.

Indications

- Provides support and immobilization in case of rotator cuff repairs, capsule joint repair, Bankart repairs or general shoulder instabilities.

Forearm length	Size	Code No.
34 - 38 cm	Small	11-0066-S
37 - 43 cm	Medium	11-0066-M
42 - 48 cm	Large	11-0066-L
+ 48 cm	X-Large	11-0066-XL

Product Composition

Polyester, Nylon, PU Foam, Plastic





Features and Benefits

Quick and effective elbow immobilisation with the option to allow or restrict motion in a certain range.

- One size-fits-most:
 - Chose LEFT or RIGHT.
 - "One push button" for telescopic length adjustment from 37 to 46 cm.
 - Eliminates holding multiple SKU's and the cost of storing multiple sizes.
- Quick adjustment of ROM:
 - By simply press and rotate settings on hinge dial.
 - Prevent patient from manipulating ROM with "lock caps" attached on press buttons. ⁽¹⁾
 - Option to lock extra in 4 positions "0° - 15° - 30° - 45°".
- Adjustable straps:
 - Convenient Y-tab closure hook for easy fitting and fast removal.

Indications

- Olecranon or distal humeral fractures (stable or internally fixed)
- Conservative treatment of elbow dislocation or subluxation
- Post-operative immobilization of elbow or ROM control for the elbow
- Severe Elbow Sprains



⁽¹⁾ ROM settings / Lock Cap



Flexibel wing cuff

	Universal	Brace length	Code No.
Left	Universal	37 - 46 cm	12-0060-1
Right	Universal	37 - 46 cm	12-0060-2

Product Composition

Polyester, Nylon, PU foam, Plastic, Fiber Glass

Elbow Strap



Features and Benefits

- The Elbow Strap is suitable for most sports activities and for daily use.
- A silicone pad provides a uniform pressure on the forearm and thus relieves the tendons.
- Strap maintains position at maximum flexion of the elbow.

Indications

- Treatment of tennis elbow (epicondylitis lateralis) and/or golfer's elbow (epicondylitis medialis)

Elbow circumference	Size	Code No.
21 - 33 cm	Universal	11-0063-U

Product Composition

Nylon, Polyester, PU Foam





D-Ring Wrist Brace & LONG D-Ring Wrist Brace

Features and Benefits

- Spacer fabric ensures a better skin aeration, providing great patient comfort.
- Anatomical design and shape specifically developed for an ideal support.
- Optimal wrist stabilization guaranteed by dual support system.
 - One volar, preformed and adjustable aluminum support.
 - Two flexible plastic stays for extra dorsal stabilisation.
- Easy application by the sleeve concept.
- Left and right version.

Indications

- Sprains of the wrist
- Tendonitis or tenovaginitis
- Overuse injuries
- Secondary stabilisation treatment after fractures
- Carpal Tunnel



	Wrist circumference	Size	Short	Long
Left	13 - 16 cm	Small	11-0043-SL	11-0045-SL
	16 - 18 cm	Medium	11-0043-ML	11-0045-ML
	18 - 20 cm	Large	11-0043-LL	11-0045-LL
	20 - 23 cm	X-Large	11-0043-XLL	11-0045-XLL
Right	13 - 16 cm	Small	11-0043-SR	11-0045-SR
	16 - 18 cm	Medium	11-0043-MR	11-0045-MR
	18 - 20 cm	Large	11-0043-LR	11-0045-LR
	20 - 23 cm	X-Large	11-0043-XLR	11-0045-XLR

Product Composition

Nylon, Polyester, Cotton, Spandex, Plastic and Aluminium

D-Ring Wrist-Thumb Brace



Features and Benefits

- Spacer fabric ensures a better skin aeration, providing great patient comfort.
- Anatomical design and shape specifically developed for an ideal support.
- A preformed, adjustable aluminium support for the wrist.
- A short and long preformed plastic support for the thumb: option for support gradation.⁽¹⁾
- Two dorsal flexible plastic stays for extra stabilisation of the wrist.
- Easy application by the sleeve concept.
- Left and right version.

Indications

- Thumb ligament injuries
- Tendonitis or tenovaginitis
- Saddle joint arthritis
- Secondary stabilisation treatment after fracture



	Wrist circumference	Size	Code No.
Left	13 - 16 cm	Small	11-0047-SL
	16 - 18 cm	Medium	11-0047-ML
	18 - 20 cm	Large	11-0047-LL
	20 - 23 cm	X-Large	11-0047-XLL
Right	13 - 16 cm	Small	11-0047-SR
	16 - 18 cm	Medium	11-0047-MR
	18 - 20 cm	Large	11-0047-LR
	20 - 23 cm	X-Large	11-0047-XLR

Product Composition

Nylon, Polyester, Cotton, Spandex, Plastic and Aluminium



Thumb Spica Splint

Features and Benefits

- Spacer fabric ensures a better skin aeration, providing great patient comfort.
- Design and shape specifically developed for an ideal support and good hand functionality.
- A symmetrical plastic insert provides a perfectly universal left-right fit.
- A compact wrist cuff doesn't obstruct your wrist mobility.
- Only 3 sizes for a reduced inventory.

Indications

- Sprain of the MCP joint (i.e. ski thumb)
- Tendonitis or tenovaginitis
- Saddle joint arthritis

	Wrist Circumference	Size	Code No.
Universal	13 - 16 cm	Small	11-0042-S
	16 - 19 cm	Medium	11-0042-M
	19 - 23 cm	Large	11-0042-L



Product Composition

Nylon, Polyester, PU Foam, Spandex, Plastic



Stack Finger Splints

Features and Benefits

- Protects the fingertip and nailbed.
- The perforated plastic allows air circulation and prevents moisture and odour buildup.
- Many sizes available.

Indications

- Injury to the extensor tendon at the DIP finger joint, also known as mallet finger (stretched or partially torn)
- After surgery of avulsion fracture of the dorsal extensor tendon



Size	Transparent	Beige
1	11-0046-1T	11-0046-1B
2	11-0046-2T	11-0046-2B
3	11-0046-3T	11-0046-3B
4	11-0046-4T	11-0046-4B
5	11-0046-5T	11-0046-5B
5 1/2	11-0046-55T	11-0046-55B
6	11-0046-6T	11-0046-6B
7	11-0046-7T	11-0046-7B

Product Composition
PVC

LOWER
EXTREMITY

*Functional
and elegant!*





Functional 4-Point Knee Orthosis

Features and Benefits

- Ultra strong and malleable aluminium structure (4 point principle).
- Contains a poly-centric hinge that imitates the anatomical movement of the knee.
- Excellent anti-slip padding prevents excessive migration of the brace.
- Numbered straps ensure a correct fitting of the brace by each patient, time after time.
NOTE: Strap 6 should be removed for the treatment of ACL injuries.
- A full set of knee condylar pads in different thickness and shapes allow good adjustment of the brace according knee width.
- Several flexion "0° - 45° - 60° - 75° - 90°" and extension stops "0° - 10° - 20° - 30° - 40°" provide the option to control the knee's range of motion adequately.

Indications

- Support after rupture or damage of ACL, PCL or MCL ligaments
- Complex instabilities of the knee joint due to trauma or degenerative reasons
- Post-Op therapy after cruciate ligament surgery

	Thigh circumference 15 cm above patella	Size	Code No.
Left	39 - 47 cm	Small	12-0022-SL
	48 - 52 cm	Medium	12-0022-ML
	53 - 60 cm	Large	12-0022-LL
	61 - 67 cm	X-Large	12-0022-XLL
Right	39 - 47 cm	Small	12-0022-SR
	48 - 52 cm	Medium	12-0022-MR
	53 - 60 cm	Large	12-0022-LR
	61 - 67 cm	X-Large	12-0022-XLR



Anti-slip padding

Product Composition

Polyester, Nylon,
PU foam, Plastic,
Fiber Glass



Universal Post-Op ROM Knee Orthosis



Features and Benefits

- A unique and versatile, universal orthosis covers a large variety of leg lengths and widths.
 - With its extendable aluminium bars and cuffs, the brace can be made extremely short (45 cm) or long (60 cm).
 - Easy circumference adjustment from quite small to XX-Large.
- After first installation remove "pull straps" so that straps attach to hook parts.
 - Hook pieces on the plastic cuff, hold the circumferential straps in place.
 - The stays remain in their position for easy re-application.
- Foam pads with anti-slip surface hold the brace well in place.
- Improved skin aeration due to limited skin coverage by the foam padding.
- Convenient click closure system facilitate removal and re-application without having to reset for each application.
- Additional foam pads ⁽³⁾ can be applied for comfort and stability at the backside of the soft thigh straps.
- Knee ROM easy to adjust or lock with dual system.
 - Determine ROM restriction by changing position of press buttons on hinge dial.
 - Prevent patient from manipulating ROM with "lock caps" ⁽¹⁾ attached on press buttons.
 - Option to lock extra in 4 positions ⁽²⁾ "0° - 15° - 30° - 45°".



Indications

- Post-operative knee immobilisation
- Stable fractures in the knee region
- Knee trauma (sprains)

	Patient height	Brace length	Code No.
Universal	1,40 - 1,90 m	44 - 60 cm	12-0021-U

Product Composition

Polyester, Nylon, PU foam, Plastic, Fiber Glass



⁽¹⁾ ROM settings / Lock Cap



⁽²⁾ Position lock



⁽³⁾ Strap foam pads



COOL ROM Knee Orthosis

Features and Benefits

- Aluminium bars on either side to provide a solid fixation.
- Hinge easy to adjust for controlled knee joint movement.
- Soft padding on the skin side for additional patient comfort.
- Convenient hook & loop closure for easy removal and application.

Indications

- Post-operative knee immobilisation
- Knee trauma (sprains)
- Minor fractures in the knee region



	Patient height	Brace length	Code No.
Universal	1,50 - 1,90 m	57 cm	12-0023-1
Short	< 1,50 m	50 cm	12-0023-2

Product Composition

Polyester, Nylon, PU Foam, Plastic, Aluminium, Low Carbon Steel

Tri-Panel Knee Immobilizer



Features and Benefits

- Two straight lateral and two contoured posterior aluminium stays.
- Stays can be removed for product cleaning and for easier re-shaping.
- Inside smooth black cotton terry towel fabric for comfort.
- Easy hook and loop strap closure for quick adjustment.
- Reduced inventory with universal sizes (thigh circumference from 35cm to 65 cm) and only 4 lengths.
- Two clever strips for keeping the main panel closed during application.

Extra for Click Closure Version

- Convenient Click Closure system for quick and easy application and removal 2 horizontal.
- 2 Horizontal silicon strips inside the immobiliser to prevent slippage.

Indications

- Post-operative knee immobilisation
- Knee trauma (sprains)
- Minor fractures in the knee region



	Patient height	Brace length	0°	20°	0° Click Closure
Universal	< 1,40 m	40 cm	11-0021-40	11-0026-40	11-0031-40
	1,40 - 1,70 m	50 cm	11-0021-50	11-0026-50	11-0031-50
	1,40 - 1,75 m	55 cm	11-0021-55	11-0026-55	11-0031-55
	> 1,70 m	60 cm	11-0021-60	11-0026-60	11-0031-60

Product Composition

Polyester, Cotton, Nylon,
PU foam, Aluminium, POM



Knee Support with Hinge

Features and Benefits

- Soft and breathable zahoprene provides comfort, and doesn't migrate.
- Elastic straps above and below the knee provide additional fixation and avoid circulation constriction.
- Extra patella support pad.
- Poly-axial hinge follows the knee movement.
- Breathable thin fabric at the knee pit provides aeration to the skin and prevents impingement of the skin.

Indications

- Knee trauma (sprains) of the collateral ligaments
- Light instabilities of the knee



Thigh circumference 15 cm above patella	Size	Code No.
36 - 39 cm	X-Small	11-0022-XS
39 - 43 cm	Small	11-0022-S
43 - 48 cm	Medium	11-0022-M
48 - 53 cm	Large	11-0022-L
53 - 58 cm	X-Large	11-0022-XL
58 - 64 cm	XX-Large	11-0022-XXL

Product Composition

Zahoprene, Polyester, Nylon, Aluminium, POM

Knee Support with Stays



Features and Benefits

- Soft and flexible zahoprene provides comfort, and doesn't migrate.
- Elastic straps above and below the knee provide additional fixation and avoid circulation constriction.
- Extra patella support pad.
- Breathable thin fabric at the knee pit provides aeration to the skin and prevents impingement of the skin.

Indications

- Moderate knee sprains
- Patella instability



Thigh circumference 15 cm above patella	Size	Code No.
34 - 39 cm	Small	11-0023-S
40 - 45 cm	Medium	11-0023-M
46 - 51 cm	Large	11-0023-L
52 - 57 cm	X-Large	11-0023-XL
58 - 64 cm	XX-Large	11-0023-XXL

Product Composition

Zahoprene, Polyester, Nylon, Aluminium, POM



Post Trauma Knee Dressing (PTKD)

The PTKD is a dressing that provides compressive support to the knee. A posterior stay gives additional immobilisation to the joint.

Features and Benefits

- Soft, comfortable PU foam, laminated with terry towel fabric provides gentle compression to the knee.
- Hook closure system allows easy application and removal, i.e. for quick wound inspection.
- The 15-mm thick foam provides already knee movement restrictions. A posterior aluminium bar makes the immobilisation even more firm
- Limited stock with just one universal size

Indications

- Intra-articular knee injuries
- Post-arthroscopic treatment of collateral and cruciate ligament injuries

Type	Brace length	Code No.
Single Pack (Plastic bag)	38 cm	11-0025-1
Six Pack (Box)	38 cm	11-0025-6

Product Composition
Polyester, PU foam,
Nylon, Aluminium



Tri-Step Ankle Lace-Up Brace



Features and Benefits

- The BraceID® Tri-Step Ankle Lace-Up is a lightweight, extremely soft and comfortable bandage.
- **Kit of stays:** more rigid anatomically shaped stays ⁽¹⁾ and straight flexible stays ⁽²⁾ offer the patient a choice between stability and comfort.
- Thin and strong figure-8 straps, in combination with a circular elastic fixation around the ankle, provide an exceptional support.
- Some spacer/3D textile elements ensure a good ventilation for the skin.
- Seamless under the heel for better comfort.
- Suits left and right feet.

Indications

- Moderate to severe ankle sprains (functional treatment)
- Prevention of ankle sprain after trauma
- Stabilisation of chronically weak ankle joints



Ankle circumference	Size	Code No.
25,5 - 28 cm	X-Small	11-0008-XS
28 - 30,5 cm	Small	11-0008-S
30,5 - 33 cm	Medium	11-0008-M
33 - 35,5 cm	Large	11-0008-L
35,5 - 38 cm	X-Large	11-0008-XL

Product Composition

Nylon, Polyester, Spandex, PP



Universal Ankle Stirrup Brace

Features and Benefits

- The BraceID® Ankle KIT is equipped with inflatable and foam cushions so that therapist or patient can choose for what best fits the circumstances.
- Air cushions are covered with a flannel textile on the skin side, which improves comfort when in contact with the skin.
- Both paddings can be easily switched by attaching them on MicroHook strips on the inside of the shells.
- Air cushions are more ideal for use in the acute phase, whereas foam pads, with a firmer density and lower profile, provide more stability for the functional phase of the rehab.

Indications

- Severe and moderate acute ankle sprains (functional and post-operative treatment)
- Stable malleolar fracture
- After care for fractures in the rehabilitative phase

Notes

- KIT = Ankle brace with all components (shells + air + foam padding)
- Easy to adjust in width to stabilise foot /ankle
- Fits easily in a classic style shoe

Description	Adult/Senior > 1,60m	Junior < 1,60m
Ankle Stirrup - AIR Padding	11-0003-UA	11-0003-JA
Ankle Stirrup - Foam Padding	11-0003-UF	11-0003-JF
Ankle Stirrup - KIT (air + foam)	11-0003-UK	Not available

Product Composition

Nylon, Polyester, EVA
& LDEP foam, POM



Ankle Wrap



Features and Benefits

- Thin bandage that fits in most shoes.
- The breathable zahoprene material ensures good ventilation of the skin.
- Two crossed elastic straps provide good support to the ankle joint. The foot insert is elastic so that the blood circulation doesn't get restricted.
- The open heel ensures sufficient freedom of movement during daily and sports activities.

Indications

- This bandage supports fragile, painful or injured ankle joints, and offers good comfort.
- Helps to prevent new injuries and is suitable for daily activities, recreation and sport.



Ankle circumference	Size	Code No.
25 - 41 cm	Universal	11-0007-U

Product Composition

Nylon, Polyester, Spandex



Shield Walkers Hi Top

Features and Benefits

- Rigid and lightweight Walker with front and back panel for more complete protection and support.
- Rocker bottom design promotes natural gait.
- Inflatable AIR liner for better fixation of foot and ankle with increased comfort. Intermittent pneumatic compression not only promotes callus formation over fracture site but also reduces local swelling.
- Supplied with additional pieces of foam padding to increase patient comfort.
- Heel wedges ⁽¹⁾ applied inside the Walker allow to position the ankle in plantar flexion (5 pieces).



Indications

- Acute ankle sprain and soft tissue injuries of the lower leg
- Stress fractures and stable or internally fixed fractures of the foot and ankle
- Achilles tendon repair

Size	Shoe size	Standard model	Air model
Small	36 - 39	12-0006-S	12-0002-S
Medium	39 - 43	12-0006-M	12-0002-M
Large	43 - 46	12-0006-L	12-0002-L
X-Large	+ 46	12-0006-XL	12-0002-XL

Product Composition

Polyester, PU foam, Nylon, Fibre Glass, Plastic



LOW Shield Walker



Features and Benefits

- Rigid and lightweight Walker with front and back panel for more complete protection and support.
- Rocker bottom design promotes natural gait.
- Inflatable AIR liner for better fixation of foot and ankle with increased comfort.
- Additional pieces of foam padding to increase patient comfort.
- Intermittent pneumatic compression not only promotes callus formation over fracture site, but also reduces local swelling.

Indications

- Immobilisation after foot or toe chirurgie
- Severe ankle sprains
- Metatarsal fractures
- Front or mid feet injuries



Size	Shoe size	Standard model	Air model
Small	36 - 39	12-0008-S	12-0010-S
Medium	39 - 43	12-0008-M	12-0010-M
Large	43 - 46	12-0008-L	12-0010-L
X-Large	+ 46	12-0008-XL	12-0010-XL



Product Composition

Polyester, PU foam, Nylon, Fibre Glass, Plastic





Walker Hi Regular

Features and Benefits

- Rigid and lightweight Walker.
- Good unrolling of the sole with intermedial plane for stance.
- Simple hook-loop closures with D-rings for tighter fixation.
- Hook protection strips at the upright bars for easier application of the liner. ⁽¹⁾
- Additional pieces of foam padding to increase patient comfort.
- Also available in SHORT version.

Indications

- Acute ankle sprain and soft tissue injuries of the lower leg
- Stress fractures and stable or internally fixed fractures of the foot and ankle
- Achilles tendon repair

Standard



Size	Shoe size	Standard model	Air model
Small	36 - 39	12-0001-S	12-0009-S
Medium	39 - 43	12-0001-M	12-0009-M
Large	43 - 46	12-0001-L	12-0009-L
X-Large	+ 46	12-0001-XL	12-0009-XL

Product Composition

Polyester, PU foam, Nylon, Fibre Glass, Plastic



Short Regular Walker



Features and Benefits

- Rigid and lightweight SHORT Walker.
- Good unrolling of the sole with intermedial plane for stance.
- Simple hook-loop closures with D-rings for tighter fixation.
- Hook protection strips at the upright bars for easier application of the liner.
- Additional pieces of foam padding to increase patient comfort.

Indications

- Immobilisation after foot or toe surgery
- Severe ankle sprains
- Metatarsal fractures
- Front or mid foot injuries



Size	Shoe size	Standard model
Small	36 - 39	12-0016-S
Medium	39 - 43	12-0016-M
Large	43 - 46	12-0016-L
X-Large	+ 46	12-0016-XL

Product Composition

Polyester, Nylon, Spandex, PU Foam, Plastic



ROM Walkers (Standard liner/AIR liner)

Features and Benefits

- Rigid and lightweight Walker with hinge for controlled movement of the ankle joint, or fixation of the ankle in a certain angle.
- Good unrolling of the sole with intermedial plane for stance.
- Inflatable AIR liner for better fixation of foot and ankle with increased comfort. The intermittent pneumatic compression not only promotes callus formation over fracture site but also reduces local swelling.
- Hook protection strips at the upright bars allow easier application of the liner.
- Additional pieces of foam padding are available to increase patient comfort.



Indications

- Acute ankle sprain and soft tissue injuries of the lower leg
- Stress fractures and stable or internally fixed fractures of the foot and ankle
- Achilles tendon repair
- Early rehabilitation of a variety of post-surgical and post-trauma indications requiring fixed or variable angular control

Size	Shoe size	Standard model	Air model
Small	36 - 39	12-0003-S	12-0004-S
Medium	39 - 43	12-0003-M	12-0004-M
Large	43 - 46	12-0003-L	12-0004-L

Product Composition

Polyester, PU foam, Nylon, Fibre Glass, Plastic





Features and Benefits

- Lightweight and solid Walker for toddlers and kids.
- Easy to apply foam padding.
- Good anchorage of the straps in the upright stays.
- Anterior and posterior shells provide a comprehensive protection and support.
- Tight fixation through the D-Ring closure.
- Efficient rocker bottom sole for an effective foot rolling.

Indications

- Greenstick fractures
- Severe sprains with skin damage (i.e. foot in wheel of a bicycle)
- Achilles Tendon release
- M. Tibialis Anterior transfer
- Toe Tip Walker – gradual correction with the aid of wedges at the forefoot


Size	Shoe size	Standard model
Small	21 - 23	12-0007-S
Medium	24 - 26	12-0007-M
Large	26 - 30	12-0007-L
X-Large	31 - 33	12-0007-XL

Product Composition

Polyester, PU foam, Nylon, Fibre Glass, Plastic



A series of horizontal dotted lines for handwriting practice, spanning the width of the page.



Find us on
social media



braceid.official



braceid.official
#braceid

Scan me!



ORTHOBROKER

Orthobroker bvba
Ter Stratenweg 11
BE-2520 Ranst
Belgium - headquarter

+32 3385 1075
info@orthobroker.com
www.braceid.com



Find us on social media 



braceid.official



braceid.official
#braceid

Wrist Embrace



	Wrist Circumference	Size	Code No.
Left	13 - 16 cm	Small	11-0040-SL
	16 - 18 cm	Medium	11-0040-ML
	18 - 20 cm	Large	11-0040-LL
	20 - 23 cm	X-Large	11-0040-XLL
Right	13 - 16 cm	Small	11-0040-SR
	16 - 18 cm	Medium	11-0040-MR
	18 - 20 cm	Large	11-0040-LR
	20 - 23 cm	X-Large	11-0040-XLR

	Wrist Circumference	Size	Code No.
Left	13 - 16 cm	Small	11-0041-SL
	16 - 18 cm	Medium	11-0041-ML
	18 - 20 cm	Large	11-0041-LL
	20 - 23 cm	X-Large	11-0041-XLL
Right	13 - 16 cm	Small	11-0041-SR
	16 - 18 cm	Medium	11-0041-MR
	18 - 20 cm	Large	11-0041-LR
	20 - 23 cm	X-Large	11-0041-XLR

Wrist-Thumb Embrace



Scan me!



ORTHOBROKER

Orthobroker bvba
Ter Stratenweg 11
BE-2520 Ranst
Belgium - headquarter

+32 3385 1075
www.braceid.com
info@orthobroker.com

Wrist Embrace

Wrist Embrace has been designed to guarantee optimum functional use and increased user comfort. By using innovative materials and design techniques, annoying or potentially painful pressure points are avoided.



Indications

- Wrist sprains
- Tendonitis / Tenovaginitis
- Overuse injuries
- Carpal Tunnel Syndrome
- Secondary treatment after fractures
- Rheumatic diseases
- Prevention
- Chronic instability of the wrist or thumb

Specific for thumb

- Rhizarthrosis
- The Quervain (Tendonitis at the thumb)



Wrist-Thumb Embrace

The Wrist-Thumb **Embrace** is a variant of the Wrist **Embrace** with an additional thumb piece. Functional and comfort aspects determined aesthetics and design decisions. A slightly more palmar positioned, ergonomic thumb splint and the "Embrace" strap technology ensure a good thumb immobilisation.



Adjustable straps

Strap-stop system



No painful pressure point on the radius




Embrace Strap

The **Embrace** strap is a new, flexible way to efficiently immobilise the thumb and is aimed at functional and comfortable rehabilitation. The **Embrace** strap system consists of two separate, adjustable straps that elegantly enclose the thumb.




Embrace strap



Strap 1 is stabilising the metacarpophalangeal joint, while **strap 2** focuses on the immobilisation of the proximal phalanx. The way in which the straps are placed around the thumb determine the degree of immobilisation.





BraceID[®]

Expert

Let's have a look !

